

SensiLance™ Pressure & Button Activated Safety Lancets

For Painless Sampling

SensiLance[™] Pressure Activated Lancets

- For capillary blood samples
- Ready to use no loading needed
- Needle automatically retracts
- Not made with natural rubber latex

SensiLance™ Button Activated Lancets

- Tri-bevel tip for virtually painless sampling
- Needle automatically retracts
- For capillary blood samples
- No lancing device required
- Not made with natural rubber latex

100	Control of the Contro	CONTRACT OF THE PROPERTY OF TH	d Color			Sterile Lancets Strait use only
	THE PROPERTY OF THE PROPERTY O	PA	Profundidad	1.8mm	2.4mm	Teer to Open Paga abrit, rasqua aqui SensiLance TM Button Activated Button Activated Button Activated
Samuel Control of Cont	Sterile Lancets Single use only - Lates Free Tear to Open Para abrir, rasgue aqui SensiLance Sensil		ección Tipo	Aguja de 28 g Aguja de 26 g	Aguja de 21 g	Button Activative Safety Lances Safety Lances In the benefit of the virtually plantists arranging Needle secondary brood samples For continuiny bood samples For continuiny bood samples For continuiny brood samples Assembly device required Assembly device required Assembly device required
	Pressure Activated Safety Lancets • For capillary blood samples • Lancet is ready to use no loading needed • Needle retracts	dynarex	Tabla de sel Pedido de reposición n.º	7125	7122	ter NO. ⁷
	automatically	21g needle 2.2mm Depth setting	^a posició	n n.º	7124	Wurning unof prosects cap Do not one support processor has been processory reached
N.	Warning: Needle is sterile only if cap is in place and lancet is undamaged. Do not re-use lancets. Lancets should be discarded in appropriate sharps containers.	Reorder No. 7112				

Reorder No.	Description	Quantity
7112	21 gauge - 2.2 mm Depth, Pressure Activated	10/100/Cs
7114	26 gauge - 1.8 mm Depth, Pressure Activated	10/100/Cs
7115	28 gauge - 1.8 mm Depth, Pressure Activated	10/100/Cs
7122	21 gauge - 2.4 mm Depth, Button Activated	10/100/Cs
7124	26 gauge - 1.8 mm Depth, Button Activated	10/100/Cs
7125	28 gauge - 1.8 mm Depth, Button Activated	10/100/Cs

7112-7125 HCPC: A4259

Other Products Available:
See our website for our
complete line of
exam gloves















dynarex[®] is your one-stop medical products brand since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com