



WE CARE LIKE FAMILY™

Sharps Containers

Dynarex high quality containers are designed to be strong and lightweight

- Containers for used medical needles and other sharp medical instruments
- Designed so needles can be dropped in the container without touching the outside
- Commonly used in EMS, hospitals, nursing homes, clinics and doctors offices, etc.
- Transparent lid allows for visual monitoring of contents to avoid overfilling

| Reorder No. | Description | Quantity |
|-------------|---|----------|
| 4622 | Sharps Container - Slide Lid, 1 qt. | 60/Cs |
| 4623 | Sharps Container - Slide Lid, 2 qt. | 60/Cs |
| 4624 | Sharps Container - Mailbox Lid, 5 qt. | 20/Cs |
| 4625 | Sharps Container - Mailbox Lid, 5.4 qt. | 20/Cs |
| 4626 | Sharps Container - Slide Lid, 1 gal. | 24/Cs |
| 4627 | Sharps Container - Slide Lid, 2 gal. | 24/Cs |
| 4628 | Sharps Container - Mailbox Lid, 3 gal. | 12/Cs |
| 4629 | Sharps Container - Slide Lid, 6 gal. | 12/Cs |
| 4630 | Sharps Shaft™, Transportable | 24/Cs |

* Depending on market demand and Dynarex stock, we will be happy to address any orders for items currently not shown.

Reorder No. 4625 - the 5.4qt container fits a BD locking bracket (#305447)



Also see our website for table paper (Reorder No. 4479 - 4482 & 4485 - 4487)

Already mobile? Scan your smart phone here for more information.



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com



WE CARE LIKE FAMILY™

Tattoo

Physician

Hospital



- Horizontal drop design for maximum capacity to fill line
- Each container locks for final disposal
- Nestable containers save space and reduce shipping costs
- Easy-to-see fill levels
- Single use containers which are disposed of with the waste inside



Reorder No. 4622

| | |
|-------------------------------------|-------|
| Sharps Container - Slide Lid, 1 qt. | 60/Cs |
| Size: 4.3"×3.5"×6.5" (11×9×16.5cm) | |



Reorder No. 4623

| | |
|-------------------------------------|-------|
| Sharps Container - Slide Lid, 2 qt. | 60/Cs |
| Size: 4.3"×3.5"×8.5" (11×9×21.5cm) | |



Reorder No. 4624

| | |
|---|-------|
| Sharps Container - Mailbox Lid, 5 qt. | 20/Cs |
| Size: 10.7"×4.6"×11.8" (27×11.5×30 cm) | |
| *Reorder No. 4624 – Compatible with Kendall Bracket Box | |



Reorder No. 4625*

| | |
|---|-------|
| Sharps Container - Mailbox Lid, 5.4 qt. | 20/Cs |
| Size: 12"×4.5"×11.3" (30.5×11.5×28.5 cm) | |
| †Reorder No. 4625 - the 5.4qt container fits a BD locking bracket (#305447) | |



Reorder No. 4626

| | |
|--------------------------------------|-------|
| Sharps Container - Slide Lid, 1 gal. | 24/Cs |
| Size: 10.3"×7"×5.5" (26×18×14 cm) | |



Reorder No. 4627

| | |
|--------------------------------------|-------|
| Sharps Container - Slide Lid, 2 gal. | 24/Cs |
| Size: 11"×7.3"×14.2" (28×18.5×36 cm) | |



Reorder No. 4628

| | |
|--|-------|
| Sharps Container - Mailbox Lid, 3 gal. | 12/Cs |
| Size: 15"×7.3"×14.2" (38×18.5×36 cm) | |



Reorder No. 4629

| | |
|--------------------------------------|-------|
| Sharps Container - Slide Lid, 6 gal. | 12/Cs |
| Size: 15.2"×10"×16" (38.5×25×40 cm) | |



Reorder No. 4630

| | |
|----------------------------------|-------|
| Sharps Shaft™ - Transportable | 24/Cs |
| Size: 6.37"×1.22" (16.17×3.1 cm) | |

*Kendall is a registered trademark of Kendall Howard
†BD is a registered trademark of Becton, Dickinson and Company.®